## Parte 2

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© 2016 by Boston Scientific Corporation or its affiliates. All Rights Reserved.

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Reference 17960317

Dear Doctor,

Thank you for contacting Boston Scientific regarding the compatibility of our Implantable devices with electromagnetic fields commonly encountered in the workplace.

#### Electromagnetic Interference (EMI) in the Workplace

Most workplace environments are unlikely to produce levels of electromagnetic energy that would interfere with the operation of implantable pacemakers and defibrillators. However, some workplace environments may contain equipment capable of producing electromagnetic energy in excess of Boston Scientific's suggested field exposure limits, and therefore may require special precautions. Environments which may produce higher levels of electromagnetic energy may include (but are not limited to) industrial manufacturing and power generation facilities.

Susceptibility to EMI could result in a temporary disruption in device function. This means that the implanted device will resume its normal mode of operation if the patient/employee either moves away from the source of EMI or turns it off. Permanent device damage as a result of EMI is very unlikely to be encountered in public and occupational environments. The following implantable device behaviors could be observed if EMI is encountered:

	Potential behaviors				
	Asynchronous pacing (intra mode)	Inhibition of pacing	Ventricular tracking	Inhibition of tachy therapy	Inappropriate tachy therapy*
Defibrillators	■	■	■	■	■
Pacemakers	■	■	■	■	■
S-ICD	■	■	■	■	■

\* Only shock therapy

Any workplace equipment that is electrically powered or contains a permanent magnet should be considered for potential EMI. The risk of EMI adversely interacting with the implantable medical device is assessed by considering numerous technical parameters such as carrier frequency, field strength, modulation scheme, duty cycle, exposure duration, and separation distance between the patient/employee's device and the field generating source.

Available resources providing guidance concerning assessment of exposure specific to workers with cardiac pacemakers and cardioverter defibrillators (ICDs) can be found in the published EN 50527-2-1: 2016 and EN 50527-2-2:2018 standards, respectively. (Note: these standards are not available

Cardiac Rhythm Management  
4100 Shadock Avenue North  
St. Paul, MN 55113-5706  
Goldnet Europe SA/NV  
A Boston Scientific Company  
Gemeinsamer Standort: 5D  
D-18111 Dagerup - Belgium  
Boston Scientific Pty Ltd  
A/Unit 40/41/42/43/44/45  
PO Box 512  
Brisbane NSW 4161 Australia  
[www.bostonscientific.com](http://www.bostonscientific.com)

Defibrillator/CRT-D	
Frequency	May Interfere with Boston Scientific defibrillators
50/60 Hz (AC)	≥ 11.7 kV/m
0 Hz (Static/DC)	≥ 1 mT or ≥ 10 G <small>Interfer designates a static magnetic flux density created by permanent magnets or from DC current flowing through a cable or source may include industrial equipment utilizing strong permanent magnets and electric motors.</small>
50 Hz (AC)	≥ 160 A/m (≥ 0.2 mT or ≥ 2 G)
60 Hz (AC)	≥ 132 A/m (≥ 0.166 mT or ≥ 1.66 G)

Interfer designates a dynamic magnetic flux density generated by current flowing through a cable or conductor carrying AC current may include high voltage power lines, transformers, open electrical panels, industrial equipment such as induction heaters, and other equipment sourced by high levels of AC current. The magnetic field strength (H-field) is derived from the

AC = Alternating Current DC = Direct Current Hz = Hertz kV/m = kilovolts per meter mT = millitesla

S-ICD	
Frequency	May Interfere with Boston Scientific Subcutaneous defibrillators
50/60 Hz (AC)	≥ 6.5 kV/m
0 Hz (Static/DC)	≥ 1 mT or ≥ 10 G <small>Interfer designates a static magnetic flux density created by permanent magnets or from DC current flowing through a cable or source may include industrial equipment utilizing strong permanent magnets and electric motors.</small>
50 Hz (AC)	≥ 160 A/m (≥ 0.2 mT or ≥ 2 G)
60 Hz (AC)	≥ 132 A/m (≥ 0.166 mT or ≥ 1.66 G)

Interfer designates a dynamic magnetic flux density generated by current flowing through a cable or conductor carrying AC current may include high voltage power lines, transformers, open electrical panels, industrial equipment such as induction heaters, and other equipment sourced by high levels of AC current. The magnetic field strength (H-field) is derived from the

AC = Alternating Current DC = Direct Current Hz = Hertz kV/m = kilovolts per meter mT = millitesla

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Reference 17960317

Dear Doctor,

Thank you for contacting Boston Scientific regarding the compatibility of our Implantable devices with electromagnetic fields commonly encountered in the workplace.

**Electromagnetic Interference (EMI) in the Workplace**

Most workplace environments are unlikely to produce levels of electromagnetic energy that would interfere with the operation of implantable pacemakers and defibrillators. However, some workplace environments may contain equipment capable of producing electromagnetic energy in excess of Boston Scientific's suggested field exposure limits, and therefore may require special precautions. Environments which may produce higher levels of electromagnetic energy may include (but are not limited to) industrial manufacturing and power generation facilities.

Susceptibility to EMI could result in a temporary disruption in device function. This means that the implanted device will resume its normal mode of operation if the patient/employee either moves away from the source of EMI or turns it off. Permanent device damage as a result of EMI is very unlikely to be encountered in public and occupational environments. The following implantable device behaviors could be observed if EMI is encountered:

	Potential behaviors					
	Asynchronous pacing (intra mode)	Inhibition of pacing	Ventricular tracking	Inhibition of tachy therapy	Inappropriate tachy therapy	Deactivation of tachy therapy*
Defibrillators	■	■	■	■	■	■
Pacemakers	■	■	■	■	■	■
S-ICD	■	■	■	■	■	■

\* Only shock therapy.

Any workplace equipment that is electrically powered or contains a permanent magnet should be considered for potential EMI. The risk of EMI adversely interacting with the implantable medical device is assessed by considering numerous technical parameters such as carrier frequency, field strength, modulation scheme, duty cycle, exposure duration, and separation distance between the patient/employee's device and the field generating source.

Available resources providing guidance concerning assessment of exposure specific to workers with cardiac pacemakers and cardioverter defibrillators (ICDs) can be found in the published EN 50527-2-1: 2016 and EN 50527-2-2:2018 standards, respectively. (Note: these standards are not available)

Defibrillatori/CRT-D	
Frequency	May Interfere with Boston Scientific defibrillators
50/60 Hz (AC)	≥ 11.7 kV/m
0 Hz (Static/DC)	≥ 1 mT or ≥ 10 G <small>Interfer designates a static magnetic flux density created by permanent magnets or from DC current flowing through a cable or source may include industrial equipment utilizing strong permanent magnets and electric motors.</small>
50 Hz (AC)	≥ 160 A/m (≥ 0.2 mT or ≥ 2 G)
60 Hz (AC)	≥ 132 A/m (≥ 0.166 mT or ≥ 1.66 G)

Interfer designates a dynamic magnetic flux density generated by current flowing through a cable or conductor carrying AC current may include high voltage power lines, transformers, open electrical panels, industrial equipment such as induction heaters, and other equipment sourced by high levels of AC current. The magnetic field strength (H-field) is derived from the following:

AC = Alternating Current DC = Direct Current Hz = Hertz A/m = Amperes per meter mT = millitesla G = Gauss kV/m = kilovolts per meter

S-ICD	
Frequency	May Interfere with Boston Scientific Subcutaneous defibrillators
50/60 Hz (AC)	≥ 6.5 kV/m
0 Hz (Static/DC)	≥ 1 mT or ≥ 10 G <small>Interfer designates a static magnetic flux density created by permanent magnets or from DC current flowing through a cable or source may include industrial equipment utilizing strong permanent magnets and electric motors.</small>
50 Hz (AC)	≥ 160 A/m (≥ 0.2 mT or ≥ 2 G)
60 Hz (AC)	≥ 132 A/m (≥ 0.166 mT or ≥ 1.66 G)

Interfer designates a dynamic magnetic flux density generated by current flowing through a cable or conductor carrying AC current may include high voltage power lines, transformers, open electrical panels, industrial equipment such as induction heaters, and other equipment sourced by high levels of AC current. The magnetic field strength (H-field) is derived from the following:

AC = Alternating Current DC = Direct Current Hz = Hertz A/m = Amperes per meter mT = millitesla G = Gauss kV/m = kilovolts per meter

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**(\*) Workers with Active Implantable Medical Devices Exposed to EMF: In Vitro Test for the Risk Assessment**

Eugenio Mattei 1,\* , Giovanni Calcagnini 1, Federica Censi 1, Iole Pinto 2, Andrea Bogi 2 and Rosaria Falsaperla 3

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# Esecuzione del test «in-vivo»

Professional Welding Station with current >200 A


AC/EMF measurements at the level of welding cable

AC/EMF measurements with patient performing welding at closest distance from welding spots and cabling.

With Real time telemetry monitoring of the cardiac device activity during testing.

Device EGM providing information for potential interference for EMF peaks transients

DC field measurement in proximity of drill magnetic holder

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© 2016 by Boston Scientific Corporation o ai affiliati. All Rights Reserved.