


Instructions for the test




The correct collection of the stool sample is an indispensable condition for obtaining good quality results.

We therefore invite you to **respect these simple instructions** and to use the material you have received (kit consisting of a test tube and plastic bag).



ON THE DAY OF THE TEST YOU SHOULD:

1. collect the faeces in a clean and dry container
2. open the test tube being careful to not spill the liquid inside
3. take out the stick (Inserted in the cap) and insert it in the faeces at various points at random (3-4 points is sufficient), taking care that a small quantity of faeces remains attached to the stick
4. reinsert once the stick in the test tube
5. close tightly
6. put the test tube in the bag and keep it in the fridge until you deliver the sample, which should be done within a maximum of two days



ATTENTION! The collection of a stool sample should not be carried out during the menstrual cycle. For the test, a specific diet is not required, you can eat normally.